


Instruction Self-Training in Clinic, Hospital, Rehab Facility (Independent Practice)



Self-training allows patients to continue structured therapy independently outside scheduled sessions - for example at the bedside, on shared facility devices, on laptops/tablets in patient rooms. This helps increase training frequency and supports better rehabilitation outcomes.

1. Setup (Therapist / Administrator)

1. Log in on the device you want to use for self-training.
2. Open **Patient Administration**.
3. Click **“Therapist / institution.”**
4. In the **Therapist data** window:
 - Activate **“Use this device for self-training.”**
 - Click **“Continue.”**




2. Patient Access (Login)


- Each patient receives a **4-digit patient ID**.
- The ID is shown in the therapy screen (top right).

How patients log in:

1. Start HeadApp → login screen appears
2. Enter **patient ID**
3. Access therapy modules directly




Therapists can switch back via **“To therapist login.”**



4. Assigning Therapy Modules (optional)

1. Go to **Patient Administration**
2. Select patient → click **“Edit”**
3. Activate **“Limit access to modules”**
4. Click **“Modules”**
5. Assign programs via double-click

👉 Only selected modules will be available to the patient




The screenshot shows the 'Patient data' form. The 'Limit access to modules' checkbox is checked, and the 'Patient is allowed to change parameters' checkbox is unchecked. The 'Modules' button is visible at the bottom right.

4.1. Restricting Patient Settings (optional)

In **Patient data**:

- Disable:
“Patient is allowed to change parameters”

👉 This ensures therapy settings remain fixed



The screenshot shows the 'Patient data' form with the 'Modules' section expanded. The 'admitted modules' list includes 'Attention and Focus Pick It', 'Time and Calendar Time It', and 'Short-term memory Flip It'. The 'inaccessible modules' list includes 'Learning and memory strategy Learn It', 'Short-term memory Pair It', 'Photo album My World', and 'Speech: Word usage'. The 'Continue' button is at the bottom right.

Key Benefits:

- Independent training outside therapy sessions
- Increased training intensity and effectiveness
- Easy login via 4-digit ID
- Full therapist control over content and settings